



MaxStar[®]

Maxstar Medical



**India Biggest Manufacturer Company
For Silicone Resuscitator & Anaesthesia Product**



Company Profile :

Welcome to Maxstar Medical

Maxstar Group of Companies is a group of three professional manufacturing companies with a focus on manufacturing and exporting of wide range of medical hospital equipment's and supplies.

We are **an ISO 9001:2008 and ISO 13485:2003 certified company**, most of our products are **CE marked**, we have three state of the art manufacturing units where we produce high standard medical products, our manufacturing units complies all the international quality standards for **Medical Devices**. Because of the Stringent quality measures adopted by the company, our products got self-introduction in the export market in a big way, today we are exporting our products up to 35 countries all over the world. The company aims to serve its customers as a **ONE STOP SUPPLIER** catering to their requirements of different products with the best of quality at the competitive pricings, using its expertise of sourcing from best available sources.

Because we strongly believe that customer relationships are at the heart of every successful company, we decided that the best use of our resources were to serve our customers in a way that no one else in the industry had approached. For that reason, we set our sights and energies on becoming a distinguished supplier of healthcare products.

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Autoclave High Pressure Horizontal

MSM-1033

Autoclave High Pressure Horizontal :-

The Autoclave High Pressure Horizontal is inner chamber made of stainless steel, door wing made of s/s Outer wall made of Stainless Steel. Three-way valve for sterilization cycle Pressure adjustable between 5 psi to 20 psi with automatic pressure control switch Water level indicator Temperature gauge, pressure release valve, safety valve. It comes also with fitted Automatic low water cut off device, Timer, Bell, Pressure breaker valve, Pressure control switch.

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Autoclave High Pressure Vertical :-

The Autoclave High Pressure Vertical is inner and outer chamber made of stainless steel. Outer chamber works as steam stock jacket. It is the most useful for sterilizing the loaded material under a steam pressure of 10 psi to 20 psi (adjustable). The space between outer wall and steam jacket is filled with high grade glass wool to minimize the temperatures loss. All sterilizers are hydraulically tested up to 40 psi. The lid made of thick stainless steel radial type closing and S.S. Basket. It is supplied complete with cord and plug to work on 220/230 volts A.C. Supply. It comes with fitted with Automatic low water cut off device, Digital temperature indicator, Dial type analog temperature gauge, Mechanical timer 60 min, Mechanical timer with alarm system.

Autoclave High Pressure Vertical

MSM-1034

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Autoclave Pressure Cooker Electric :-

The Autoclave Pressure Cooker Electric is made of Aluminum / Stainless Steel & Mirror Finish. It comes with fitted Pressure Gauge, Steam Release Valve & Safety Valve.

Capacity (Approx.): 6 Ltr. to 56 Ltr..

Autoclave Pressure Cooker Electric

MSM-1036

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Autoclave Pressure Cooker

MSM-1035

Autoclave Pressure Cooker :-

The Autoclave Pressure Cooker is made of Aluminum / Stainless Steel & Mirror Finish. It comes with fitted Pressure Gauge, Steam Release Valve & Safety Valve.

Capacity (Approx.): 6 Ltr., 9 Ltr., 12 Ltr., 14 Ltr. & 16 Ltr.

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Autoclave Win Nut Type Electric :-

The Autoclave Pressure Cooker is made of Aluminum / Stainless Steel & Mirror Finish. It comes with fitted Pressure Gauge, Steam Release Valve & Safety Valve.

Capacity (Approx.): 6 Ltr., 9 Ltr., 12 Ltr., 14 Ltr. & 16 Ltr.

Autoclave Win Nut Type Electric

MSM-1037

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com

Autoclave Sterilizer



Hot Air Sterilizer

Hot Air Sterilizer (Oven):-

The Hot Air Sterilizer (Oven) is double walled Inside Chamber of Aluminum / Stainless Steel. & Outside M.S. painted with attractive stove enamel with adjustable wire mesh nichol plated trays, on/off Rocker switch with Indicator & Thermostat Indicator. It works on 220/230 Volts A.C. The Temperature Thermostatically controlled by hydraulic thermostat from 50°C to 250°C with a sensitivity of +/- 1°C in Bare Oven. Beaded heating elements are placed in ribs at the bottom and sides. Built in L-shape thermometer.

MSM-1040

info@maxstarmedical.com

www.maxstarmedical.com
www.resuscitatorindia.com



Maxstarmedical



Facebook/maxstarmedical



twitter.com/medisafedelhi



Contact Us:-

Get in Touch : Mr. Akash Pandey (Sales & Marketing Dpt.)

Address : A 10-11, Bhandari House 1st Floor Dr.Mukherjee
Nagar, Delhi-110009 (INDIA)



+(91) - (11)-27650687
+(91) - 9810309116



info@maxstarmedical.com
sales@maxstarmedical.com



www.maxstarmedical.com