



MaxStar[®]

Maxstar Medical



**India Biggest Manufacturer Company
For Silicone Resuscitator & Anaesthesia Product**



Company Profile :

Welcome to Maxstar Medical

Maxstar Group of Companies is a group of three professional manufacturing companies with a focus on manufacturing and exporting of wide range of medical hospital equipment's and supplies.

We are **an ISO 9001:2008 and ISO 13485:2003 certified company**, most of our products are **CE marked**, we have three state of the art manufacturing units where we produce high standard medical products, our manufacturing units complies all the international quality standards for **Medical Devices**. Because of the Stringent quality measures adopted by the company, our products got self-introduction in the export market in a big way, today we are exporting our products up to 35 countries all over the world. The company aims to serve its customers as a **ONE STOP SUPPLIER** catering to their requirements of different products with the best of quality at the competitive pricings, using its expertise of sourcing from best available sources.

Because we strongly believe that customer relationships are at the heart of every successful company, we decided that the best use of our resources were to serve our customers in a way that no one else in the industry had approached. For that reason, we set our sights and energies on becoming a distinguished supplier of healthcare products.

 **Air Cushions**



MSM-1058

Air Cushions :-

Indian Rubber Air Rings are being exported all over the world in huge quantities.

Square Air Rings can also be made against bulk requirement.

Please send us your list of requirements and we shall be glad to send our C.I.F. rates and samples.

BP Machine Aneroid



MSM-1043

BP Machine Aneroid

A compact and reliable Sphygmomanometer supplied with adult Velcro fabric cuff and precision manometer with clip for fix.

Features:-

- Emulsoid Bladder & Soft emulsoid Bulb.
- Soft Cotton / Nylon cuff.
- Leather bag easy to carry.
- Latex-free for safety.
- Certified Accurate to +/- 3mmHg.
- Laser-engraved dial face for accuracy.

BP Machine Mercury



MSM-1041

BP Machine Mercury

There is a button beside the LCD screen, if you put the button down, it will be Oscillometric measuring method. It's like digital B.P. meter, it's used by puff the bulb manually. When the LCD column reach to about 160-180mmHg, stop to puff, about 10seconds, the value will display on the screen. It's easier than Stethoscopy method.

LCD Mercury Free Sphygmomanometer is a new generation mercury-free blood-pressure meter which is still made with Krotkoff Sound.

BP Machine Stand



MSM-1042

BP Machine Stand

Mercurial Sphygmomanometer is traditional blood pressure measuring instruments with mercury locking facility, Range 0-300 mm Hg packed in lightweight metal case.

Features:- Easy to read with clear tube
Accurate reading.
Emulsoid Bladder & Soft emulsoid Bulb
Mercury lock facility.

The Height of pole is adjustable (with range of 105mc to 146cm) Five wheels with brakes

 **BP Machine Wall**



MSM-1136

BP Machine Wall



Breast Pump



MSM-1044

Breast Pump

To understand mother's needs & baby's behavior, we closely work with leading universities in the field of breastfeeding and human lactation. Helping mothers to express milk is our core competency.

Manual breast pumps are by far the most affordable and portable pump options available. If you only plan to express milk occasionally then a manual pump may be the best choice for you. A manual pump gives you complete control over the speed and suction of the pump.

You can easily take a manual breast pump with you anywhere since it does not require electricity, is compact, lightweight and easily fits in a diaper bag or purse.

 **Corrgated Sheet**



Corrgated Sheet

MSM-1140

E.N.T Set



E.N.T Set

BAYONET MOUNTED LOCKING HEAD SYSTEM, XBH 2.5 V ILLUMINATION CHROME PLATED AND KNURLED HANDLE, REUSABLE SPECULA 1, 2 & 3 RHEOSTAT FOR ADJUSTMENT OF ILLUMINATION.

Ent Complete Diagnostic Set With Ophthalmoscope
1 Light weight C size Handle with Rheostat Button
1 Otoscope Head with 3 Reusable Ear Speculums
1 Ophthalmoscopes Head with Selection Wheel and Contain 24 Optical Lens
1 Bent Arm Illuminator, 1 Nasal Speculum
1 Tongue Depressor / Blade, 1 Tongue Depressor Blade Holder,
2 Laryngeal mirrors

MSM-1045

 **ECG Blubs**



ECG Blubs

We offer high quality ECG Blubs. Rubber bulb fitted 30 mm in diameter fitted with electrode cup. Bulb colors available in Red, Black and Light Blue.

MSM-1059

ECG Machine 6 Channel



ECG Machine 6 Channel

MSM-1148

ECG Machine



MSM-1046

ECG Machine

ECG is a transthoracic interpretation of the electrical activity of the heart over time captured and externally recorded by skin electrodes. It is a noninvasive recording produced by an ECG device.

Features:

- Digital isolation technology and signal processing solution, digital filter
- 7" TFT, touch screen, led backlight
- Patient management, name and age and id.
- Detail analysis report.

 **Enema Syringe**



Enema Syringe

MSM-1170



Foetal Doppler



Foetal Doppler

The Foetal Doppler is a hand-held ultrasound detector used to detect the heartbeat of a fetus for prenatal care. A Foetal Doppler is called a baby Doppler, is used to detect Fetal Heart Rate (FHR) after about 12 weeks pregnancy. It is used in hospitals, clinics or at home.

- Easy to read LCD with on/off back light.
- Clear digital fetal heart rate display.
- Batteries uses size 2 AA size batteries.

MSM-1064

 **Forceps**



MSM-1060

Forceps

Specifications:-Length: 0.7/1.2/1.6/1.8/2.3mm

Jaw size: 1.8mm/ 2.3mm /3.0mm.

Coating: With /without.

Spike: With /without needle.

Specifications:-O. D. 3.0mm *1600/1800/2300mm.

O. D. 2.3mm *1600/1800/2300mm.

O. D. 1.8mm *700/1200/1600/1800/2300mm.

Hot Water Bottles



Hot Water Bottles

- Fill warm water up to two third of the capacity.
- Unscrew the stopper.
- Expel the air by pressing the bottle gently.
- Put the stopper and tight it.
- Bottle is ready for use.
- Keep one third of the bottle empty.
- Do not use boiling water.
- After filling press gently to expel air.
- Stopper should not be loose.
- Avoid mishandling by children.
- When not in use keep the empty bottle up side down after removng the stopper.

MSM-1047

 **ICE Bags**



MSM-1048

ICE Bags

We offer high quality & long life ICE Bags. A convenient reusable double usage Gel pack suitable for musculoskeletal pain and high fever conditions.

Infusion Pump



MSM-1062

Infusion Pump

Features:

- Compact structure
- Designed with microcomputer technology
- Has multi-alarm function and takes on the function of speedy transportation simultaneously
- Installs rechargeable battery can be use for rescue work in field, in verge areas and satisfied with all kind of needs of hospital.
- Mainly use as phleboclysis and chemotherapy infusion in the department of internal medicine, surgery, emergency, pediatrics, gynecology and department of blood, ICU monitor ward and CCU monitor ward.

 **Maekintosh Sheeting**

Maekintosh Sheeting



MSM-1171



Otoscope



MSM-1050

Otoscope

Our Otoscope has LED Light (White Light) with high focus. Our Otoscope makes use of a scratch resistant variable focus, optical quality GLASS Lens with a 10X magnification. It produces Bright White Light. The LED lighting has longlife meaning you will do notneed to replace a bulb for long time. Our Otoscope LED'S produce a white light like sunshine and do not have typical 'yellow' spectrum light produced by halogen bulbs for more natural color rendition of skin lesions. Our LED lighting also uses very minimal battery power. Our Otoscope is built from Chrome Plated Brass which allows for a beautiful finish along with longlife.



Oxygen Concentrator



MSM-1063

Oxygen Concentrator

It has a Base Substructure and Recessed Casters which protects it against damage and breakage.

It is equipped with a Rotary Poppet Valve, delivering up to 5 lpm of therapeutic oxygen in home, hospital, and clinical settings. With hundreds of thousands of units in daily use throughout the world, this low power consumption concentrator is known for high performance, easy maintenance, and unmatched reliability. (Flow lockdown and other unique features make this a concentrator of choice among today's patients and equipment providers.)

 **Percussion Hammers**

Percussion Hammers



MSM-1173



Pulse Oximeter



MSM-1061

Pulse Oximeter

We offer Pulse Oximeter is used to rate heart beat of persons. Pulse oximeters are medical devices that monitor the level of oxygen in a patient's blood and alert the health-care worker if oxygen levels drop below safe levels, allowing rapid intervention. These devices are essential in any setting in which a patient's blood oxygen levels requires monitoring like operations, emergency and intensive care, and treatment and recovery in hospital wards. It is used in Hospitals, Clinic & Nursing Homes.



Ring Pessary Silicon & Rubber



MSM-1072

Ring Pessary Silicon & Rubber

The Ring Pessary is made of high grade raw material. It is 100% Latex Free & Autoclavable at 134°C for 15 minutes (As per International Standard).

Available Size:-50 mm, 55 mm, 60 mm, 65 mm, 70 mm, 75 mm, 80 mm, 85 mm, 90 mm, 95 mm, 100 mm, 105 mm, 110 mm.

Available Color:-Red, Blue, White (As per customer's requirement)

 **Rubber Blulbs**



Rubber Blulbs

MSM-1174

 **Rubber Tubes**



Rubber Tubes

MSM-1051

 **Stethoscopes Tubes**



MSM-1057

Stethoscopes Tubes

The Stethoscope Tubing is made of high grade raw material. The Stethoscope Tubing is available in many colors.

 **Stethoscopes**



Stethoscopes

We offer good quality stethoscope, Reliable & Clear acoustics offer an advantage to Health Care Professionals when listening to Heart, Lung and Blood Pressure Sounds, Light Weight· Chrome Plated Chest Piece, Brass / Aluminum / S. Steel, Excellent Acoustic Diaphragm, External Non- Chill Rim, Y Shaped tubing in Variety of colors, Soft sealing Ear tips.

MSM-1175

Syringe Distrage



Syringe Distrage

MSM-1176

 **Syringe Pump**



Syringe Pump

MSM-1177



Thermometer Digital



Thermometer Digital

The Digital Thermometer is accurate and fast Thermometer which is environmental friendly and very safe to use. Clear vision of temperature and beep alert sound further add to its features. It can be used orally, rectally & under armpit

Then, the temperature taken last time would be display with a "M" mars in the right corner.

When L and flashing °C (°F) display, the thermometer is ready for temperature measuring. But if the room temperature is higher than 32.0°C (90°F), the room temperature would be displayed instead of "L" The beep will sound when temperature measuring is complete. Press On/OFF button of the unit. (The unit will automatically power off after approximately eight minute after use.

MSM-1178

Thermometer Mercury



Thermometer Mercury

The Thermometer Mercury is accurate. It can be used orally, rectally & under armpit.

Specification:-

Clinical Economy Model Oral,

Prismatic type Lens front.

Dual scale Cylindrical bulb Range: 35-45°C & 95-108°F

Accuracy $\pm 0.10^{\circ}\text{C}$ & $\pm 0.2^{\circ}\text{F}$

Packed in plastic case.

MSM-1054



Tuning Forks



MSM-1055

Tuning Forks

The tuning fork is made of high grade stainless steel. It has been used for a wide variety of experimental and exhibitory purposes. These were used as frequency standards in research and lecture demonstrations.

 **Vaginal Speculums**



Vaginal Speculums :-

MSM-1056



Contact Us:-

Get in Touch : Mr. Akash Pandey (Sales & Marketing Dpt.)

Address : A 10-11, Bhandari House 1st Floor Dr.Mukherjee
Nagar, Delhi-110009 (INDIA)



+(91) - (11)-27650687
+(91) - 9810309116



info@maxstarmedical.com
sales@maxstarmedical.com



www.maxstarmedical.com